









# IN ADVANCE

Melt YPD in microwave and pour plates.

Restreak "Golden Yeast" onto YPD to isolate single white colonies\*\*

# DAY OF LAB

1. Label 1 microfuge tubes: **NO DNA**.
2. Add 500  $\mu$ l of **EZ - Solution I** to the **NO DNA** tube.
3. Use a sterile pipet tip, toothpick or inoculating loop, scrape an isolated white colony of yeast off the petri dish.
4. Swirl the colony into the **EZ - Solution I** in the **NO DNA** tube.
5. Discard the pipet tip, toothpick or inoculating loop into a waste receptacle to be decontaminated.
6. Microfuge the tube you have prepared, coordinating with another group or setting up a "blank" microfuge tube with 500  $\mu$ l of water. Spin the tubes for 30 seconds at full speed.
7. After microfuging the tubes, the cells will have collected as a white pellet in the bottom of the tube. Remove as much of the supernatant as you can, using a pipet and discarding the liquid into a waste receptacle to be decontaminated.
8. Resuspend the pellet in 150  $\mu$ l of **EZ - Solution II**, pipetting up and down to make a homogeneous solution.
9. Label two more microfuge tubes: **Control DNA**, **crtYB\* DNA**. Add 5  $\mu$ l of the appropriate DNA to the tubes, changing tips between aliquots. Skip this step if the DNA has already been aliquoted for you.



**10.** Add 50  $\mu$ l of cells from the **NO DNA** tube to the **Control DNA** tube. Pipet up and down to mix.

**11.** Add 50  $\mu$ l of cells from the **NO DNA** tube to the **crtYB\* DNA** tube. Pipet up and down to mix.

**12.** Add 500  $\mu$ l of **EZ - Solution III** to each microfuge tube. Solution III will be goopy, but the amount you pipet does not need to be precise. Pipet up and down to mix, changing tips between tubes.

**13.** Incubate the tubes at 30° Celcius (C) for one hour. Periodically flick the tubes to mix during the incubation.

**14.** Label the media side of three SC-trp petri dishes as **NO DNA/white colony**, **Control DNA/white colony**, and **crtYB\* DNA/white colony**. Add your initials and today's date to each.

**15.** Pipet 200  $\mu$ l of each sample onto the media of the appropriate petri dish. Spread evenly across the dish with a sterile spreader. \*\* Discard spreader and microfuge tubes into the waste receptacle to be decontaminated.

**16.** Incubate petri dishes, media side up, for 2 days at 30° C.

After the petri dishes have incubated for 2 days, count the colonies of each color in every dish.

